



Anti-Spna2 polyclonal antibody (CABT-BL3430)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to residues surrounding amino acids 162 of mouse Spna2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Chicken, Human, Mouse, Rat
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Buffer	In PBS, pH 7.2 (30% glycerol, 0.5 mg/mL BSA, 0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Entrez Gene ID	20740
Protein Refseq	NP_001171139.1
UniProt ID	B2RXX6