



Anti-SPNS3 (C-terminal) polyclonal antibody (CABT-BL3431)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from the C terminal region of mouse SPNS3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane.
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	77577
Protein Refseq	NP_084208

UniProt ID

[Q9D232](#)
