



Anti-SPPA polyclonal antibody (CABT-BL3433)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant sppA protease IV from Arabidopsis thaliana
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Arabidopsis thaliana, C. reinhardtii
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cell inner membrane; Multipass membrane protein
Format	Liquid
Buffer	Whole serum
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Encodes a putative protease SppA (SppA).
--------------	--

GENE INFORMATION

Entrez Gene ID	843737
----------------	------------------------

Protein Refseq [NP_565077](#)

UniProt ID [Q9C9C0](#)
