



Anti-SRSF12 (aa 179-228) polyclonal antibody (CABT-BL3441)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to a region within internal sequence amino acids 179-228 (RSRSKSLQKR SKSIGKSQSS SPQKQTSSGT KSRSHGRHSD SIARSPCKSP) of Human SRRP35.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [135295](#)

Protein Refseq [NP_542781](#)

UniProt ID [Q8WXF0](#)
