



Anti-SSNA1 (full length) polyclonal antibody (CABT-BL3445)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length protein, corresponding to amino acids 1-119 of Human SSNA1
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, ICC/IF
Cellular Localization	Nuclear. Centrosome.
Format	Liquid
Buffer	PBS, pH 7.2
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	8636
Protein Refseq	NP_003722

UniProt ID

[A0A024R8G6](#)
