



# Anti-SST polyclonal antibody (CABT-BL3448)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Somatostatin (6-28)
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	N/A
<b>Purification</b>	IgG fraction
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Dot, RIA, WB
<b>Cellular Localization</b>	Secreted
<b>Format</b>	Liquid
<b>Buffer</b>	Liquid antiserum
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

<b>Introduction</b>	Somatostatin (also known as growth hormone-inhibiting hormone (GHIH) or somatotropin release-inhibiting factor (SRIF)) or somatotropin release-inhibiting hormone[citation needed] is a peptide hormone that regulates the endocrine system and affects neurotransmission and cell proliferation via interaction with G protein-coupled somatostatin receptors and inhibition of the release of numerous secondary hormones. Somatostatin regulates insulin and glucagon.
---------------------	---

# GENE INFORMATION

Entrez Gene ID [24797](#)

---

Protein Refseq [NP\\_036791](#)

---

UniProt ID [P60042](#)

---