



Anti-SST polyclonal antibody (CABT-BL3449)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to amino acids 1-14 of cyclic Human Somatostatin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	IHC-P
Cellular Localization	Secreted
Format	Liquid
Buffer	1% BSA, Whole Serum, 0.01M PBS
Preservative	0.1% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction	Somatostatin (also known as growth hormone-inhibiting hormone (GHIH) or somatotropin release-inhibiting factor (SRIF)) or somatotropin release-inhibiting hormone[citation needed] is a peptide hormone that regulates the endocrine system and affects neurotransmission and cell proliferation via interaction with G protein-coupled somatostatin receptors and inhibition of the release of numerous secondary hormones. Somatostatin regulates insulin and glucagon.
---------------------	---

GENE INFORMATION

Entrez Gene ID	6750
Protein Refseq	NP_001039
UniProt ID	P61278