



Anti-STAMBPL1 (aa 1-50) polyclonal antibody (CABT-BL3455)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 1-50 (MDQPFTVNSLKKLAAMPDHT DVSLSPEERV RALSKLGCNI TISEDITPRR) of Human STAMBPL1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	AMSH-like protease is an enzyme that in humans is encoded by the STAMBPL1 gene.
--------------	---

GENE INFORMATION

Entrez Gene ID [57559](#)

Protein Refseq [NP_065850](#)

UniProt ID [Q96FJ0](#)
