



Anti-Stk30 (C-terminal) polyclonal antibody (CABT-BL3458)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant GST fusion protein corresponding to C-terminus 111 amino acids of mouse Stk30.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Buffer	In PBS (0.1% proclin, 2% Block Ace)
Preservative	None
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Entrez Gene ID	26448
Protein Refseq	NP_036103.1
UniProt ID	Q9WVS4