



Anti-STX8 (aa 1-236) polyclonal antibody (CABT-BL3465)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length protein, corresponding to amino acids 1-236 of Human Syntaxin 8.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cell Membrane; single-pass type IV membrane protein
Format	Liquid
Size	100 µg
Buffer	1X PBS, pH 7.2
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction The gene is a member of the syntaxin family. The encoded protein is involved in protein

trafficking from early to late endosomes via vesicle fusion and exocytosis. A related pseudogene has been identified on chromosome 12. Alternative splicing results in multiple transcript variants.

GENE INFORMATION

Entrez Gene ID	9482
----------------	----------------------

Protein Refseq	NP_004844
----------------	---------------------------

UniProt ID	Q9UNK0
------------	------------------------
