



Anti-STX8 (cytoplasmic domain) polyclonal antibody (CABT-BL3466)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant cytoplasmic domain of rat syntaxin 8 (aa2- 215).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	WB, IP, ICC
Preservative	None
Storage	-20°C, avoid repeated freeze/thaw cycles

BACKGROUND

Introduction	The gene is a member of the syntaxin family. The encoded protein is involved in protein trafficking from early to late endosomes via vesicle fusion and exocytosis. A related pseudogene has been identified on chromosome 12. Alternative splicing results in multiple transcript variants.
--------------	--

GENE INFORMATION

Entrez Gene ID	9482
----------------	----------------------

Protein Refseq [NP_004844](#)

UniProt ID [Q9UNK0](#)
