



Anti-STXBP5L polyclonal antibody (CABT-BL3468)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide as a part of mouse Stxbp5l conjugated to an immunogenic carrier protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	IgG fraction
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-Fr
Reconstitution	Reconstitute in 500 µl of sterile water. Centrifuge to remove any insoluble material.
Cellular Localization	Cytoplasm (Potential). Cell membrane; Peripheral membrane protein (Potential). Membrane; Peripheral membrane protein (Potential). Note=Cytoplasmic, and associated with vesicular membranes and the plasma membrane (Potential).
Format	Lyophilised
Buffer	PBS, Whole serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [207227](#)

Protein Refseq [NP_001108083](#)

UniProt ID [Q5DQR4](#)
