



Anti-SULT6B1 (aa 180-229) polyclonal antibody (CABT-BL3474)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: AEFLGFFLTG EQIQTISVQS TFQAMRAKSQ DTHGAVGPFL FRKGEVGDWK, corresponding to C terminal amino acids 180-229 of Human SULT6B1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasmic
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Store at -20°C. Stable for 12 months at -20°C

GENE INFORMATION

Entrez Gene ID	391365
Protein Refseq	NP_001027549

UniProt ID

[Q6IMI4](#)
