



Anti-SUSD3 (aa 36-85) polyclonal antibody (CABT-BL3475)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide from within a region of N terminal amino acids 36-85 (RLPPQATFQV LRGNGASVGT VLMFRCPSNH QMVGSGLLTC TWKGSIAEWS) of Human SUSD3
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Reconstitution	Add 50 µl of distilled water for a final concentration of 1 mg/ml.
Cellular Localization	Membrane; Single pass membrane protein
Format	Lyophilised
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [203328](#)

Protein Refseq [NP_659443](#)

UniProt ID [Q96L08](#)
