



# Anti-SWI5 (aa 109-138) polyclonal antibody (CABT-BL3478)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	KLH conjugated synthetic peptide designed between amino acids 109-138 from the central region of Human SWI5
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Cellular Localization</b>	Nucleus.
<b>Format</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at 4°C (up to 6 months). For long term storage store at -20°C

## GENE INFORMATION

<b>Entrez Gene ID</b>	<a href="#">375757</a>
<b>Protein Refseq</b>	<a href="#">NP_001035100</a>

UniProt ID

[Q1ZZU3](#)

---