



# Anti-SYNGR4 (N-terminal) polyclonal antibody (CABT-BL3484)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide: HIPKSLQELANSC, corresponding to amino acids 2-14 of Human Synaptogyrin 4
<b>Isotype</b>	IgG
<b>Source/Host</b>	Goat
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, WB
<b>Cellular Localization</b>	Membrane.
<b>Format</b>	Liquid
<b>Buffer</b>	0.5% BSA, 5mg/ml Tris, pH 7.3
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

<b>Introduction</b>	This gene encodes an integral membrane protein. The gene belongs to the synaptogyrin gene family. Like other members of the family the protein contains four transmembrane regions. The exact function of this protein is unclear. [provided by RefSeq, Jul 2008]
---------------------	---

# GENE INFORMATION

Entrez Gene ID [23546](#)

---

Protein Refseq [NP\\_036583](#)

---

UniProt ID [O95473](#)

---