



# Anti-SYNGR4 (N-terminal) polyclonal antibody (CABT-BL3485)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	A 14 amino acid peptide from near the amino terminus of Human Synaptogyrin 4 (CAA09754).
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Mouse, Rat, Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, IHC-P
<b>Cellular Localization</b>	Membrane.
<b>Format</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at 4°C, stable for one year. Should not be exposed to prolonged high temperatures. Avoid freeze / thaw cycles.

## BACKGROUND

<b>Introduction</b>	This gene encodes an integral membrane protein. The gene belongs to the synaptogyrin gene family. Like other members of the family the protein contains four transmembrane regions. The exact function of this protein is unclear. [provided by RefSeq, Jul 2008]
---------------------	---

# GENE INFORMATION

Entrez Gene ID [23546](#)

---

Protein Refseq [NP\\_036583](#)

---

UniProt ID [O95473](#)

---