



Anti-SYS1 (N-terminal) polyclonal antibody (CABT-BL3490)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to N-terminal residues of human SYS1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Golgi apparatus membrane; Multi-pass membrane protein.
Format	Liquid
Buffer	50% Glycerol, PBS
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	SYS1 forms a complex with ADP-ribosylation factor-related protein ARFRP1 (MIM 604699) and targets ARFRP1 to the Golgi;apparatus (Behnia et al., 2004 (PubMed 15077113)).
---------------------	--

GENE INFORMATION

Entrez Gene ID [90196](#)

Protein Refseq [NP_001184058](#)

UniProt ID [Q8N2H4](#)
