



Anti-SYTL1 (aa 123-171) polyclonal antibody (CABT-BL3491)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 123-171 (EAEAAVKEKE EGPEPRLTID EAPQERLRET EASDP EEASQ AQEDPGQGDQ) of Human SYTL1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cell membrane. Peripheral membrane protein tightly bound to the cytoplasmic side of cellular membranes.
Format	Liquid
Buffer	PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction Synaptotagmin-like protein 1 is a protein that in humans is encoded by the SYTL1 gene.

GENE INFORMATION

Entrez Gene ID [84958](#)

Protein Refseq [NP_001180237](#)

UniProt ID [Q8IYJ3](#)
