



Anti-TADA1 (aa 288-335) polyclonal antibody (CABT-BL3501)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 288-335 (REVIPTHTVY ALNIERIITK LWHPNHEELQ QDKVHRQRLA AKEGLLLC) of Human TADA1L.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	TADA1L is a protein subunit of the human STAGA complex (SPT3; (MIM 602947)/TAF9 (MIM 600822)/GCN5 (MIM 602301);acetyltransferase complex), which is a chromatin-modifying multiprotein complex (Martinez et al., 2001 (PubMed;11564863)).
---------------------	---

GENE INFORMATION

Entrez Gene ID [117143](#)

Protein Refseq [NP_444281](#)

UniProt ID [Q96BN2](#)
