



Anti-TAF10B polyclonal antibody (CABT-BL3502)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fusion protein (Fruit fly (<i>Drosophila melanogaster</i>))
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Fruit fly
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	33468
Protein Refseq	NP_477418
UniProt ID	Q9XZT7