



Anti-TAF3 (aa 540-589) polyclonal antibody (CABT-BL3503)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region from internal sequence amino acids 540-589 (RDREREKDKNKDKSKEKDKVKEKEKDKETGRETYPWKE FLKEEEADPYK) of Human TAF3
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	The highly conserved RNA polymerase II transcription factor TFIID (see TAF1; MIM 313650) comprises the TATA box-binding;protein (TBP; MIM 600075) and a set of TBP-associated
---------------------	---

factors (TAFs), including TAF3. TAFs contribute to promoter;recognition and selectivity and act as antiapoptotic factors (Gangloff et al., 2001 (PubMed 11438666)).

GENE INFORMATION

Entrez Gene ID	83860
----------------	-----------------------

Protein Refseq	NP_114129
----------------	---------------------------

UniProt ID	Q5VWG9
------------	------------------------
