



# Anti-TAF4B polyclonal antibody (CABT-BL3508)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                    |   |
|--------------------|---|
| Immunogen          | Recombinant protein corresponding to amino acids of human TAF4B.                                    |
| Isotype            | IgG   |
| Source/Host        | Rabbit  |
| Species Reactivity | Human   |
| Purification       | Antigen affinity purification   |
| Conjugate          | Unconjugated  |
| Applications       | IHC, IF   |
| Format             | Liquid  |
| Buffer             | In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)   |
| Preservative       | 0.02% Sodium Azide  |
| Storage            | Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing. |

## BACKGROUND

|              |  |
|--------------|--|
| Introduction | TATA binding protein (TBP) and TBP-associated factors (TAFs) participate in the formation of the TFIID protein complex, which is involved in initiation of transcription of genes by RNA polymerase II. This gene encodes a cell type-specific TAF that may be responsible for mediating transcription by a subset of activators in B cells. |
|--------------|--|

# GENE INFORMATION

Entrez Gene ID [6875](#)

---

Protein Refseq [NP\\_005631](#)

---

UniProt ID [J3KTH2](#)

---