



Anti-TAS2R104 (cytoplasmic domain) polyclonal antibody (CABT-BL3514)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide from a cytoplasmic domain of mouse TAS2R4, conjugated to an immunogenic carrier protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-Fr
Cellular Localization	Membrane; Multi-pass membrane protein.
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [387340](#)

Protein Refseq [NP_996894](#)

UniProt ID [Q7M723](#)
