



## Anti-TAS2R108 (cytoplasmic domain) polyclonal antibody (CABT-BL3515)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide from a cytoplasmic domain ( C-terminal) of Rat TAS2R4 conjugated to an immunogenic carrier protein.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Rat
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, IHC-P, IHC-Fr
<b>Cellular Localization</b>	Membrane; Multi-pass membrane protein.
<b>Format</b>	Liquid
<b>Buffer</b>	Whole serum
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles.

### GENE INFORMATION

Entrez Gene ID [554302](#)

Protein Refseq [NP\\_001019857](#)

---

UniProt ID [Q67ET0](#)

---