



Anti-TAS2R19 polyclonal antibody (CABT-BL3517)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide as a part of human TAS2R48 conjugated to immunogenic carrier protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-Fr
Cellular Localization	Membrane.
Format	Liquid
Buffer	PBS, Whole serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Taste receptor type 2 member 19 is a protein that in humans is encoded by the TAS2R19 gene.
--------------	---

GENE INFORMATION

Entrez Gene ID	259294
Protein Refseq	NP_795369
UniProt ID	P59542