



Anti-TAS2R19 (internal region) polyclonal antibody (CABT-BL3518)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthesized peptide derived from an internal region of Human TAS2R48.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane.
Format	Liquid
Buffer	PBS, 0.88% Sodium chloride, 50% Glycerol pH: 7.40
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

BACKGROUND

Introduction Taste receptor type 2 member 19 is a protein that in humans is encoded by the TAS2R19 gene.

GENE INFORMATION

Entrez Gene ID [259294](#)

Protein Refseq [NP_795369](#)

UniProt ID [P59542](#)
