



# Anti-TAS2R30 polyclonal antibody (CABT-BL3519)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	A synthetic peptide as a part of human TAS2R47 conjugated to an immunogenic carrier protein
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, IHC-P, IHC-Fr
<b>Cellular Localization</b>	Membrane; Multipass membrane protein.
<b>Format</b>	Liquid
<b>Buffer</b>	PBS, Whole serum
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C short term (1-2 weeks). Add glycerol to a final volume of 50%, aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction** Taste receptor type 2 member 30 is a protein that in humans is encoded by the TAS2R30 gene.

# GENE INFORMATION

Entrez Gene ID [259293](#)

---

Protein Refseq [NP\\_001091112](#)

---

UniProt ID [P59541](#)

---