



Anti-TBC1D3 (C-terminal) polyclonal antibody (CABT-BL3529)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	This antibody is specific for the C Terminus Region of the target protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Buffer	20 mM Potassium Phosphate, 150 mM Sodium Chloride, pH 7.0
Preservative	None
Storage	Store at -20°C. Avoid freeze-thaw cycles.

BACKGROUND

Introduction	Acts as a GTPase activating protein for RAB5. Does not act on RAB4 or RAB11.
--------------	--

GENE INFORMATION

Entrez Gene ID	729873
----------------	------------------------

Protein Refseq [NP_001116863.2](#)

UniProt ID [A0A087X195](#)
