



Anti-TBC1D9 (C-terminal) polyclonal antibody (CABT-BL3534)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region between residue 1118 and the C-terminus (residue 1168) of human TBC1D9.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µl
Buffer	8mM PBS, 60mM Citrate, 150mM Tris, pH 7-8
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction	TBC1 domain family member 9 is a protein that in humans is encoded by the TBC1D9 gene.
--------------	--

GENE INFORMATION

Entrez Gene ID [23158](#)

Protein Refseq [NP_055945](#)

UniProt ID [Q6ZT07](#)
