



Anti-TCF15 (aa 2-51) polyclonal antibody (CABT-BL3543)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: AFALLRPVGA HVLYPDVRL SEDEENRSES DASDQSFGCC EGPEAARRGP, corresponding to amino acids 2-51 of Human TCF15
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, WB
Reconstitution	Add 50 µl of distilled water for a final concentration of 2mg/ml.
Cellular Localization	Nuclear
Format	Lyophilised
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction The protein encoded by this gene is found in the nucleus and may be involved in the early

transcriptional regulation of patterning of the mesoderm. The encoded basic helix-loop-helix protein requires dimerization with another basic helix-loop-helix protein for efficient DNA binding.

GENE INFORMATION

Entrez Gene ID	6939
----------------	----------------------

Protein Refseq	NP_004600
----------------	---------------------------

UniProt ID	A4LBB6
------------	------------------------
