



Anti-TCHP (aa 306-370) polyclonal antibody (CABT-BL3546)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, from a region within internal sequence amino acids 306-370 of Human TCHP (Q9BT92)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Cellular Localization	Cytoplasmic, cytoskeleton
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Trichoplein keratin filament-binding protein is a protein that in humans is encoded by the TCHP gene.
---------------------	---

GENE INFORMATION

Entrez Gene ID [84260](#)

Protein Refseq [NP_001137324](#)

UniProt ID [A0A024RBM9](#)
