



Anti-TCOF1 (aa 4-193) polyclonal Antibody (CABT-BL3548)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to N terminal amino acids 4-193 of Human Treacher Collins syndrome protein
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus; nucleolus.
Format	Liquid
Buffer	PBS
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	May be involved in nucleolar-cytoplasmic transport. May play a fundamental role in early embryonic development, particularly in development of the craniofacial complex.
---------------------	--

GENE INFORMATION

Entrez Gene ID	6949
Protein Refseq	NP_000347
UniProt ID	Q13428