



Anti-TDH (internal region) polyclonal antibody (CABT-BL3551)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from an internal domain of mouse TDH
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	IHC-P, IHC-Fr, WB
Cellular Localization	Mitochondrial
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	58865
Protein Refseq	NP_067455

UniProt ID

[Q8K3F7](#)
