



Anti-TDRD6 (aa 1700-1800) polyclonal antibody (CABT-BL3552)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH derived from within residues 1700 - 1800 of Human TDRD6.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasmic; localises to chromatoid bodies in spermatids
Format	Liquid
Size	100 µg
Buffer	1% BSA, PBS, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [221400](#)

Protein Refseq [NP_001010870](#)

UniProt ID [O60522](#)
