



Anti-TECTB (aa 61-91) polyclonal antibody (CABT-BL3554)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH conjugated synthetic peptide designed between amino acids 61-91 from the N-terminal region of Human TECTB
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cell membrane. Secreted; extracellular space; extracellular matrix. Found in the non-collagenous matrix of the tectorial membrane.
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

BACKGROUND

Introduction	The genes for alpha-tectorin (MIM 602574) and beta-tectorin encode the major noncollagenous proteins of the tectorial membrane of the cochlea.
---------------------	--

GENE INFORMATION

Entrez Gene ID	6975
Protein Refseq	NP_478129
UniProt ID	Q96PL2