



Anti-TEX37 polyclonal antibody (CABT-BL3556)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant_protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Purified
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Lyophilized Powder
Buffer	PBS with 0.02% Sodium Azide, PH7.2
Preservative	0.02% Sodium Azide
Storage	-20°C, Avoid freeze / thaw cycles

BACKGROUND

Introduction	TEX37 (Testis Expressed 37) is a Protein Coding gene.
--------------	---

GENE INFORMATION

Entrez Gene ID	74221
Protein Refseq	NP_083101.1
UniProt ID	Q9DAG4