



Anti-TH polyclonal antibody (CABT-BL3566)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length SDS denatured protein (purified from pheochromocytoma) (Rat).
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Rat
Purification	Protein A purified
Conjugate	Unconjugated
Applications	IHC-Fr, IHC-P, IHC-P, IHC-FoFr, WB, IP, ICC
Format	Liquid
Size	100 µl
Buffer	50% Glycerol, PBS, 100µg/ml BSA, 150mM Sodium chloride, 10mM HEPES, pH 7.5
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	catalyzes the rate limiting step in catecholamine biosynthesis; may be involved in response to hypoxia [RGD, Feb 2006]
---------------------	--

GENE INFORMATION

Entrez Gene ID [7054](#)

Protein Refseq [NP_000351](#)

UniProt ID [P07101](#)
