



Anti-Th polyclonal antibody (CABT-BL3570)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Native purified rat Th.
Source/Host	Sheep
Species Reactivity	Rat
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB, IHC, IF
Format	Liquid
Buffer	In 10 mM HEPES, 260mM NaCl, pH 7.5 (50% glycerol, 10% BSA)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	catalyzes the rate limiting step in catecholamine biosynthesis; may be involved in response to hypoxia [RGD, Feb 2006]
--------------	--

GENE INFORMATION

Entrez Gene ID	25085
Protein Refseq	NP_036872

UniProt ID

[P04177](#)
