



Anti-TH (aa 65-255) polyclonal antibody (CABT-BL3572)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Purified recombinant protein of human tyrosine hydroxylase (aa 65 - 255).
Source/Host	Pig
Species Reactivity	Human
Conjugate	Unconjugated
Applications	IHC
Preservative	None
Storage	-20°C, Avoid freeze / thaw cycles

BACKGROUND

Introduction	catalyzes the rate limiting step in catecholamine biosynthesis; may be involved in response to hypoxia [RGD, Feb 2006]
--------------	--

GENE INFORMATION

Entrez Gene ID	7054
Protein Refseq	NP_000351
UniProt ID	P07101