



Anti-THOC4 polyclonal antibody (CABT-BL3575)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region between residue 75 and 125 of Human Aly (AAH52302.1).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IP
Cellular Localization	Cytoplasmic and Nuclear
Format	Liquid
Buffer	Tris citrate/phosphate buffer, pH 7-8
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction	The protein encoded by this gene is a heat stable, nuclear protein and functions as a molecular chaperone. It is thought to regulate dimerization, DNA binding, and transcriptional activity of basic region-leucine zipper (bZIP) proteins.
---------------------	--

GENE INFORMATION

Entrez Gene ID [10189](#)

Protein Refseq [NP_005773](#)

UniProt ID [E9PB61](#)
