



# Anti-TIGD6 polyclonal antibody (CABT-BL3586)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Recombinant protein corresponding to amino acids of human TIGD6.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Antigen affinity purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, IHC
<b>Format</b>	Liquid
<b>Buffer</b>	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Introduction</b>	The protein encoded by this gene belongs to the tigger subfamily of the pogo superfamily of DNA-mediated transposons in humans. These proteins are related to DNA transposons found in fungi and nematodes, and more distantly to the Tc1 and mariner transposases. They are also very similar to the major mammalian centromere protein B.
---------------------	---

## GENE INFORMATION

Entrez Gene ID [81789](#)

---

Protein Refseq [NP\\_001230182](#)

---

UniProt ID [Q17RP2](#)

---