



Anti-TLX1NB (aa 72-121) polyclonal antibody (CABT-BL3591)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide to a region within: PHPLNKQKHH PPHPSQTQKD LVPRSPQLEK SRIRLRRTLRL NLGGGRGQRG, corresponding to C terminal amino acids 72-121 of Human TDI.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Entrez Gene ID	100038246
Protein Refseq	NP_001078867