



Anti-TM4SF20 (aa 107-156) polyclonal antibody (CABT-BL3592)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 107-156 (QALLKGPLMC NSPSNSNANC EFSLKNISDI HPESFNLQWF FNDSCAPPTG) of Human TM4SF20.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cell membrane.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [79853](#)

Protein Refseq [NP_079071](#)

UniProt ID [Q53R12](#)
