



Anti-TMC2 polyclonal antibody (CABT-BL3595)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide from mouse TMC2, conjugated to an immunogenic carrier protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-Fr
Cellular Localization	Cell Membrane; Multi-pass membrane protein
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Transmembrane channel-like protein 2 is a protein that in humans is encoded by the TMC2 gene.
---------------------	---

GENE INFORMATION

Entrez Gene ID [192140](#)

Protein Refseq [NP_619596](#)

UniProt ID [Q8R4P4](#)
