



Anti-TMEM145 (internal region) polyclonal antibody (CABT-BL3606)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from the internal region of human TMEM145
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Membrane; Multi-pass membrane protein.
Format	Liquid
Buffer	50% Glycerol, PBS (without Mg ²⁺ and Ca ²⁺), 150mM Sodium chloride, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

GENE INFORMATION

Entrez Gene ID	284339
Protein Refseq	NP_775904

UniProt ID

[Q8NBT3](#)
