



Anti-TMEM146 (aa 36-85) polyclonal antibody (CABT-BL3607)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|------------------------------|---|
| Immunogen | Synthetic peptide within residues: LIQDVQGDRL YFHPTTTRLI KHPCEKNIAL YLGKQVFFTM DNFETSLLPF, corresponding to N terminal amino acids 36-85 of Human TMEM146 |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Human |
| Purification | Immunogen affinity purified |
| Conjugate | Unconjugated |
| Applications | WB, ELISA |
| Cellular Localization | Membrane; Single-pass type I membrane protein |
| Format | Liquid |
| Buffer | 2% Sucrose, PBS |
| Preservative | None |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. |

GENE INFORMATION

| | |
|-----------------------|---------------------------|
| Entrez Gene ID | 257062 |
| Protein Refseq | NP_689997 |

UniProt ID

[Q86XM0](#)
