



Anti-TMEM163 polyclonal antibody (CABT-BL3610)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal region amino acids 143-192 (AAVHSAHREYIACVILGVIFLLSSICIVVKAIHDLSTRL LPEVDDFLFSV) of Human TMEM163
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	81615
Protein Refseq	NP_112185

UniProt ID

[Q8TC26](#)
