



Anti-TMEM184A (C-terminal) polyclonal antibody (CABT-BL3613)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide from the C terminus of Human TMEM184A, conjugated to KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-P
Cellular Localization	Membrane.
Format	Liquid
Buffer	99% PBS
Preservative	0.02% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction	TMEM184A (transmembrane protein 184A) is a protein-coding gene. An important paralog of this gene is TMEM184C.
---------------------	--

GENE INFORMATION

Entrez Gene ID [202915](#)

Protein Refseq [NP_001091089](#)

UniProt ID [Q6ZMB5](#)
