



Anti-TMEM2 (aa 1-50) polyclonal antibody (CABT-BL3615)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 1-50 (MYAAGSRGHS PAFLQPQNGN GHRSPGYVPG KVVPLRPAPP PKNHASAKLT) of Mouse TMEM2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane.
Format	Liquid
Buffer	97% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [83921](#)

Protein Refseq [NP_001028931](#)

UniProt ID [Q5FWI3](#)
